**Nelson** deluxe

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### B

Betriebsart 5, 15, 20  
 BICOAG – 7, 13, 15, 30  
 BICUT 1 – 7, 13, 15, 30  
 BICUT 2 15  
 BICUT 2 – 7, 13, 30  
 Bipolares Koagulieren 7, 8, 9, 10, 15  
 Bipolares Schneiden 7, 8, 9, 10, 13, 14, 15  
 CONTACTCOAG 6, 13, 15, 30  
 FORCEDCOAG 6, 13, 15, 30  
 MICROCUT 1 – 6, 13, 15, 30  
 MICROCUT 2 – 6, 13, 15, 30  
 Monopolares Koagulieren 6, 8, 10, 14, 15  
 Monopolares Schneiden 6, 8, 10, 13, 14, 15  
 POWERCUT 1 – 6, 13, 15, 30

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

**Symbols on *Nelson*<sup>deluxe</sup> (GN 640)**

	Neutral electrode
	Foot pedal
	Attention, see instructions for use
	Neutral electrode for HF with insulated earth grounding connection
	Type CF unit with defibrillation protection
	Potential equalization
	Fuse
	Non-ionizing current
	Alternating current
	Power ON/OFF switch
	Labelling of electric and electronic devices according to directive 2002/96/EC (WEEE), see Disposal

## Contents

1.	Safe handling. ....	3
1.1	Safe handling in conformance with IEC regulations .....	3
2.	Product description. ....	5
2.1	Components required to operate the device. ....	5
2.2	Intended use .....	5
2.3	Mode of operation .....	7
3.	Preparation and setup. ....	10
3.1	Stacking devices .....	10
4.	Working with the device .....	11
4.1	System set-up .....	11
4.2	Operating the device. ....	13
4.3	Functional test .....	21
5.	Cleaning, disinfecting, sterilization and product care. ....	22
5.1	Disinfecting and cleaning .....	22
5.2	Care and checking for damage .....	22
5.3	Sterilization method and parameters .....	22
6.	Maintenance .....	22
6.1	Safety inspection. ....	22
7.	Recognizing and remedying errors .....	24
7.1	Calling up the activation, error and On Phase memories ....	24
7.2	Calling up the menu Service Functions. ....	25
7.3	Error messages on the display .....	25
7.4	Changing fuses .....	29
8.	Technical Service. ....	29
9.	Disposal .....	29
10.	Accessories and spare parts .....	29
11.	Technical specifications .....	30
12.	Index. ....	30

## 1. Safe handling

These instructions for use are not suitable for instructing a beginner in surgical techniques with high frequency devices. This documentation does not contain a general description of the applications of high frequency surgical devices (see the O.R. manual and other relevant literature).

➤ Only allow persons to use the product and its accessories who possess the requisite training, knowledge and experience to do so.

- Prior to operating the device, check thoroughly to ensure that it is working properly and that it is not damaged in any way.
- To prevent damage caused by improper assembly or use and in order not to compromise the warranty and liability:
  - Use the product only according to these instructions for use.
  - Follow the safety and maintenance instructions.
  - Only combine Aesculap products with each other.
  - Use the device in conformance with prevailing standards (see Safe handling in conformance with IEC regulations).
- Accessories must be inspected regularly: Electrode cables and endoscopic accessories, in particular, must be checked for possible damage to the insulation.
- Keep the instructions for use accessible for the O.R. personnel.
- Adhere to prevailing standards.

### 1.1 Safe handling in conformance with IEC regulations

#### Operating environment

Certain applications of HF devices can cause sparks!

- Do not use the device in an area where there is a risk of explosion.
- When performing operations on the cranial and thoracic regions, either avoid the use of explosive anaesthetic gases and inflammable gases (e.g. laughing gas, oxygen) or suck away these substances.
- Only use non-flammable cleaning and disinfecting agents and solvents (for adhesives).
- If flammable cleaning and disinfecting agents or solvents are used: Be sure that these substances have evaporated before using HF surgery.
- Be sure that no inflammable liquids accumulate under the patient's body or in body cavities (e.g., the vagina). Before using the HF device, wipe up all fluids.
- Be sure there are no endogenous gases in the O.R. environment that could ignite.
- Be sure that any materials that were soaked in liquid oxygen (e.g., cotton wool, gauze) are kept far enough away from the HF environment that they cannot ignite.

#### Risk of damage to other devices!

HF devices create potentially harmful magnetic fields during normal use.

- Be sure that no electronic devices that could be influenced by an electromagnetic field are set up in the vicinity of the HF device.

## Patient safety

### Danger from inadequate preparation or operational errors in the HF device!

- Be sure that the HF device is in perfect working order.
- Be sure that no conductive fluids (e.g. blood, amniotic fluid) have entered the foot pedal or the hand switch.
- Be sure there is no short circuit in the cable of the foot pedal or the hand switch.

### Risk of the patient incurring burns arising from accidental activation of the HF device!

- Switch off the HF device immediately, if it was activated accidentally, using the OFF switch .
- Always exercise extreme caution when using the foot pedal and/or the hand switch.

Always be sure to do the following:

- Bed the patient in such a way that s/he is not in contact with any metal parts that are grounded or have a significant electric capacity against ground (e.g. operating table, fixtures). If necessary, interpose antistatic drape.
- Be sure that the patient does not come in contact with any damp cloths, drape or bedding.
- Lay anti-static drape against parts of the body that tend to develop large amounts of sweat in order to safeguard against skin-to-skin contact with the torso.
- Siphon off urine with a catheter.
- When operating on the open heart, ground the HF device via the potential equalization device.
- Prior to performing HF surgery on patients who wear pacemakers, consult with the attending cardiologist in order to safeguard against irreparable damage to, or failure of, the pacemaker device.
- If possible, remove from the vicinity of HF electrodes any electrodes from physiological monitoring devices that are not equipped with protective resistors or HF dampers.
- Do not use needle electrodes for intraoperative monitoring.
- Arrange the wires and cables of monitoring devices in such a way that they do not come in contact with the patient's skin.
- Keep the pathways of the wires and cables of high-frequency electrodes as short as possible, and arrange them in such a way that they do not touch the patient or each other.

- To avoid unwanted coagulation in other areas of the body when performing surgery on small parts of the body and on such high-resistance areas as bones and joints: Use bipolar technique in those areas.
- Always set the power output of HF devices to as low a level as possible.
- In the event that the standard settings do not yield satisfactory results, be sure that:
  - neutral electrodes are properly emplaced.
  - the operational electrodes are clean.
  - devices are plugged in and connected properly.
- Set the acoustic signal that sounds when an electrode is activated at a volume that makes the signal easy to hear.
- Never put down an active electrode on or next to the patient.
- Set aside active electrodes that are not needed in such a way that they do not come in contact with the patient.
- In operations where continuous contact between the electrodes and the patient cannot be avoided (e.g. endoscopic operations), actuate OFF switch  15 immediately to switch off the HF device after an accidental activation.
- Do not remove hot electrodes from the patient's body immediately after performing cutting or coagulation procedures.

### Risk of injury to the patient due to an unintended rise of the HF output power in case of any fault in the HF device!

- Stop using the device immediately if even the slightest functional anomaly is detected.

## Application of the neutral electrode

In the monopolar HF technique, the neutral electrode serves to conduct the current introduced to the body at the operation site back to the HF device.

Typically, either reusable neutral electrodes made of conductive rubber, or self-adhesive disposable electrodes are usually used. Both of these types of electrodes are quite suitable for use on uneven body surfaces.

To avoid a temperature rise at the exit point of the current, the following conditions must be ensured:

- Contact surfaces of adequate size between neutral electrodes and the body.
- High electrical conductivity between the neutral electrodes and the body.

To prevent burns caused by the neutral electrode, always ensure the following:

- Select areas of application for neutral electrodes that result in the shortest possible electrical pathways between active and neutral electrodes; and that are longitudinal or diagonal to the body (since muscles with the same orientation as the fibrilla are more conductive).
- In operations in the thorax area, make certain the path of the current does not run across the chest and does not run through the heart.
- Depending on the operation field, apply the neutral electrode at the nearest upper arm or thigh, if possible.  
When applying disposable self-adhesive electrodes, follow the instructions of the supplier with regard to the area of application.
- Ensure that the application area is free of both scar tissue and bony outgrowth.
- Make sure that the current pathway does not intersect any implants (e.g., bone nails, bone plates, endoprotheses).
- Prior to applying the neutral electrode:
  - Remove bushy hair growth.
  - Clean the area to be operated on (do not use alcohol, as it dries out the skin and increases contact resistance).
  - In patients with poor circulation, massage or brush the application area.
- Apply neutral electrodes systematically over the entire surface. Secure reusable neutral electrodes with rubber bands or elastic bandage, thereby ensuring that they do not come loose when the patient moves. At the same time, ensure that such measures do not constrict blood flow (risk of necroses).
- Apply neutral electrode type 3 in such a way that the black, conductive surfaces are in contact with the patient and the blue surface faces away from the body.
- To forestall raised contact resistance caused by epidermal dehydration arising from the use of damp cloths or electrode paste: do not under any circumstances use damp cloths or electrode paste.
- To prevent burns arising from short bursts of high electric current: Ensure that no fluids (e.g. irrigation fluids, disinfectants, blood, urine) get between the patient and the neutral electrode.
- To prevent fluid buildup, reduced blood flow or perfusion, pressure necroses, displacement of the neutral electrodes and contact resistance: do not lay neutral electrodes under patient's buttocks or back.
- Check to make sure that the electrical pathway of the HF device does not intersect any EEG electrodes.

## 2. Product description

### 2.1 Components required to operate the device

The following components are required to operate the high-frequency surgical device *Nelson*<sup>deluxe</sup> (GN 640):


- Power cord
- Foot pedal
- Neutral electrode
- Electrode handle with "fingertip" buttons.
- Operating electrode
- Connecting cable
- Bipolar equipment

### 2.2 Intended use

The *Nelson*<sup>deluxe</sup> (GN 640) is intended for use in surgery, for generating the HF power for monopolar and bipolar incisions and coagulating in micro- and macro-surgery.











#### Use of the monopolar and bipolar modes of operation

All operational modes are illustrated by two different pictograms, one of which is a simplified pictogram, which is always visible, and the other a detailed pictogram, which is displayed when the operational mode is selected.





Each of these pictograms displays a typical application and its effects.

## Monopolar incisions


		<b>MICROCUT 1 + 2</b> Spark-controlled incising with dynamically minimized power (stronger for MICROCUT 1; dynamic range determined by the dosage setting): Incisions virtually independent of electrode shape and tissue properties
		
		<b>POWERCUT 1</b> Powerful incisions with a low degree of scabbing: Incisions in normal or low-ohm tissues – voltage regulation
		<b>POWERCUT 2–4</b> Adjustable degree of scabbing as an additional feature compared to POWERCUT 1 – voltage regulation (maximal for POWERCUT 4)
		
		
		<b>SEALCUT</b> Incisions with alternate cutting pulses and coagulation phases (e.g. for the low-hemorrhage removal of polyps with a loop electrode, or for papillotomy)

## Monopolar coagulation

		<b>SOFTCOAG</b> Gentle, tissue-preserving coagulation, not sinking into the tissue (e.g. for critical coagulations with resectoscopes and for use with monopolar endoscopic instruments)  Additional functions auto-COAG off and Vario-function are not supported
		<b>CONTACTCOAG</b> Coagulation over broad surfaces with clamps or forceps; hemostasis of oozing hemorrhaging; cauterization of large areas of issue (e.g., in tumor surgery); coagulation of small surfaces  Additional functions supported: autoCOAG off, Vario-function
		<b>FORCEDCOAG</b> Superficial coagulation with spark formation  Additional functions auto-COAG off and Vario-function are not supported
		<b>SPRAYCOAG</b> Surface coagulation through spark discharge, no contact between tissue and electrode; special application: Hemostasis in difficult-to-access bone crevices and in parenchymatous tissues; optimized coagulation is achieved with supplemental application of the argon function  Additional functions auto-COAG off and Vario-function are not supported

## Note

When selecting the coagulation mode **CONTACTCOAG**, the additional functions supported are displayed as well.

For a detailed description of the additional function autoCOAG off see Special settings (Menu Special Functions) and for the Vario-function see Activating the high frequency.

### Bipolar techniques



In bipolar incisions and coagulation, good incision and coagulation results can be achieved with low power settings.

Since neutral electrodes need not be applied, the risk to the patient associated with them is eliminated.


#### Note

With the bipolar technique, optimum results (especially in minimally invasive surgery) can only be achieved with special instruments, see Aesculap brochure C-304-81.

### Bipolar incisions

		<b>BICUT 1</b> Incisions with bipolar instruments (e.g. Aesculap BITOM instruments)
		<b>BICUT 2</b> Incisions with bipolar scissors (e.g. Aesculap Metzenbaum scissors PM 400 R)

### Bipolar coagulation

		<b>BICOAG</b> Dosage 0.1–19.5: Microsurgical application with fine instrument tips, tissue sparing coagulation with strong depth effect Dosage from 20: Coagulation with the entire range of bipolar coagulation instruments  Additional functions supported: autoCOAG on/off; power-on delay, Vario-function
---	---	---

## 2.3 Mode of operation

The high-frequency surgical device *Nelson*<sup>deluxe</sup> (GN 640) is microprocessor-controlled and transforms the mains voltage into a high-frequency alternating current for monopolar and bipolar application.

In monopolar operation, the device can be activated either by a handle with fingertip buttons or through a pedal control; in bipolar mode, the device is activated through an instrument with fingertip buttons, through a foot pedal or via the autoCOAG function.

It is not possible to activate several outputs simultaneously.

### Data memory functions

With the push of a button, the user can create their own device settings and call up empirical data previously entered into the system (see Working with the device).

As on any computer, the data saved remains in the device memory even when the device is switched off.

### Monitoring functions

#### Self-test

When switched on, the device automatically conducts a self-test that assesses operational components, the warning signal, the microprocessor and hardware functionality.

#### Ongoing tests while in operation

Safety-relevant functions and signals are cyclically tested while the device is in operation. If any fault is detected, the HF generator switches off. The display shows an error message (see Error messages on the display). The device makes a distinction between operational errors and device errors.

### Neutral electrode monitoring

The neutral electrode monitoring function minimizes the risk of burns at the neutral electrode application site. When an error is detected, the monopolar HF generator switches off and cannot be activated.

It is possible to continue working in bipolar mode.

The following events are monitored:



- Undivided neutral electrode (type 1):  
When the connection between a neutral electrode and the HF device is severed, or the plug is disconnected, an error message appears on the display.
- Divided disposable neutral electrode (type 2)  
– or –  
Reusable divided neutral electrode with isolated outer surfaces (type 3):  
When a divided neutral electrode is used, the PDM system (Permanent Dynamic Monitoring System) also checks if the neutral electrode is in correct contact with the patient. If an error is detected (e.g., the neutral electrode is partially detached from the patient), an error message appears on the display.

## Output diagrams


Output ( $P_{out}$ ) in relation to dosage:

Operational modes of monopolar CUT and COAG




Output ( $P_{out}$ ) in relation to dosage:

Operational modes of bipolar CUT and COAG





## Adaptation curves



Operational modes of monopolar CUT and COAG  $P = f(R_L)$

Operational modes of bipolar CUT and COAG:  $P = f(R)$ ;  
(Measurement with original cable:  $C_L = 175 \text{ pF}$ )



## Maximum peak output voltage ( $U_p$ )


These graphs are intended to enable the user to judge the suitability of the HF surgical device or its output setting for a given accessory, concerning the quality of the insulation. This judgment has to be based on the peak output voltage given in the documentation accompanying the accessory.

Maximum peak output voltage ( $U_p$ ) correlated with dosage:  
Monopolar operational mode

Maximum peak output voltage ( $U_p$ ) correlated with dosage:  
Bipolar operational mode

## 3. Preparation and setup

Aesculap accepts no responsibility in the event the following parameters are not adhered to.

- When setting up and using the product, adhere to:
  - government rules and regulations regarding installation and operation,
  - and national regulations on fire and explosion protection.

### 3.1 Stacking devices

#### Note

Materials that serve to connect the high-frequency surgical device **Nelson<sup>deluxe</sup>** (GN 640) to other devices (with and without stacking devices) are available as accessories and spare parts, see *Accessories and spare parts*.

- Do not exceed a maximum stack height of 475 mm.
- Make certain the support (e.g. table, ceiling support, equipment trolley) is sufficiently stable.
- Aesculap devices without integrated stacking devices: Remove the feet of the upper device by hand or with the aid of a screwdriver.
- Place the upper device on the one below it.
- Put the stacking plates in place (see Fig.).
- Check to ensure that the devices are properly seated in one other by lifting them slightly.
- Position devices so that they are stable.

#### Moving stacks of devices

- Always lift stacks from underneath the bottommost device.

## 4. Working with the device

### 4.1 System set-up


#### Plugging the device in and switching it on

The voltage of the main power supply must correspond to the voltage indicated on the rating plate at the back of the device.

The device is outfitted with a universal power supply unit, which means that it automatically adapts to voltages ranging from 100 V to 240 V and that no switching needs to be performed.

- Plug the power cord into the power connection **21** at the back of the control unit.
- Then insert the unit plug into the room power supply socket.
- Turn on the control unit with the power ON **16** switch.  
The power ON **17** indicator lamp is illuminated. The HF device is performing a self-test.


On switching on the device for the first time after the self-test, the basic factory settings are shown on the display.

Each time the device is switched on thereafter, the device settings used before the last switch-off are displayed.

#### Using the function keys

The adjustment of all settings of the HF device **Nelson<sup>®</sup> deluxe** (GN 640), e.g. dosage setting, operational mode selected or menu settings, can be carried out via the function keys **3, 4, 5, 6** and **14**.


The actuation of a function key is confirmed by an acoustic signal when the key is pressed.

#### Connecting the neutral electrode

The connection socket neutral electrode **12** has been designed for a 6.3 mm jack plug.




When no neutral electrode is plugged in, or in the case of a neutral electrode fault, the neutral electrode symbol lights up red . If the neutral electrode has been connected correctly, the neutral electrode symbol appears green .



### Connecting active electrodes

Various plugs can be connected to the active output sockets **7, 8, 9** and **10**. The permissible plug-socket combinations are listed in the following.

#### Output sockets HF – bipolar 1 and 2

Plug types		
Coaxial plug	2-pin plug (pin spacing: 29 mm)	Plug: 4 mm
		

#### Output sockets HF – monopolar 1 and monopolar 2

Plug types		
Coaxial plug	3-pin plug	Single plug: 4 mm
		

WARNING

When plugging in the monopolar 4-mm plugs, care has to be taken that the correct socket is used; otherwise, error 48  (pin error monopolar) will be triggered. The red ring marks the socket carrying the HF current.

DANGER

Risk for patients and users due to accidental current supply to an instrument!

- Make certain that not more than 1 patient cable is connected to any output socket.


### Connecting the foot pedal

The foot pedal plug and the connection socket foot pedal **11** are marked with an arrow each. To be able to connect the plug to the connection socket correctly, these two arrows must point to each other.

### User settings (Menu User Settings)

In the User Settings menu, the following settings can be changed:


- Type of neutral electrode
  - Type 1: non-divided neutral electrode;  
Display shows 
  - Type 2: divided neutral electrode;  
Display shows 
  - Type 3: Reusable, divided neutral electrode (GN 320);  
Display shows 
- Acoustic signal volume
  - Acoustic signal during activation: Range 40–100 %
  - Warning signal: no settings option, always 100 %
- Date/Time
  - Date in the format: DD.MM.YY
  - Time in the format: hh:mm
- Maximum HF switching time
  - 30 seconds (default setting)
  - 45 seconds
  - 60 seconds
  - 20 s for automatic switch-off with autoCOAG Short (fixed value)
- Press function key  **2**.  
The display shows:

#### Note


After a waiting time of 30 seconds, or when the user presses function key , the system returns to the main menu.

- Press the function key  6.  
The menu User Settings is displayed:

- Change the setting selected by actuating function keys  14, as shown on the display.
- Confirm the setting by pressing the function key  and return to the main menu.  
- or -
- To exit the menu without adopting the changes:  
Press the function key  or wait for 30 seconds.

## 4.2 Operating the device

### Selecting and de-selecting HF outputs

The HF outputs of bipolar 1 and 2 as well as monopolar 1 and 2 modes can be selected and de-selected at will in accordance with application.

- To select and de-select HF outputs: Press the corresponding function key  3, 4, 5, 6.  
The menu card of HF output that can be operated is displayed, underlaid in gray.  
The functions of function keys  14 are allocated to the gray menu card.
- Set function parameters such as Dosage setting, operational mode and activation mode, using function keys  14.

### Selecting CUT/COAG mode

Various operational modes can be selected for each output that is activated.

Available operational modes CUT:

- Monopolar output:
  - MICRO CUT 1–2
  - POWER CUT 1–4
  - SEAL CUT
- Bipolar output:
  - BICUT 1
  - BICUT 2


Available operational modes COAG:

- Monopolar output:
  - SOFT COAG
  - CONTACT COAG
  - FORCED COAG
  - SPRAY COAG
- Bipolar output:
  - BICOAG


## Selecting operational mode CUT

- Press the function key . The menu bar operational mode CUT is displayed:

## Monopolar

## Bipolar


- Press the function key  repeatedly until the desired CUT mode is highlighted in black.
- Confirm the setting by pressing the function key  and return to the main menu.
- or -
- To exit the menu without adopting the changes: Press the function key  or wait for 30 seconds.

## Selecting operational mode COAG

- Press the function key . The menu bar operational mode COAG is displayed:

## Monopolar

- Press the function key  repeatedly until the desired COAG mode is highlighted in black.
- Confirm the setting by pressing the function key  and return to the main menu.
- or -
- To exit the menu without adopting the changes: Press the function key  or wait for 30 seconds.

## Note

If an output is allocated with an operational mode that is not supported by the selected activation mode autoCOAG off, the system automatically selects a manual method of activation.

## Selecting a dosage

A dosage must be set for each HF output that is selected.

The allowable dosage range depends upon the operational mode selected:


Operational mode		Dosage min.	Dosage max.
Monopolar	MICROCUT 1	5	300
	MICROCUT 2	5	300
	POWERCUT 1	5	300
	POWERCUT 2	5	250
	POWERCUT 3	5	200
	POWERCUT 4	5	150
	SEALCUT	5	100
	SOFTCOAG	5	100
	CONTACTCOAG	1	120
	FORCEDCOAG	1	120
	SPRAYCOAG	1	120
Bipolar	BICUT 1	1	100
	BICUT 2	1	100
	BICOAG	0,1	Foot pedal: 100 autoCOAG: 50

The size of the graduations between dosages depends upon operational mode and dosage range:

Operational mode	Dosage range	Graduation
Monopolar CUT	1–50	1
	55–100	5
	110–300	10
Monopolar COAG	1–50	1
	55–120	5
Bipolar CUT	1–50	1
	55–100	5
Bipolar COAG	0,1–10	0,1
	10,5–20	0,5
	21–50	1
	55–100	5

- To select HF output for which the dosage is to be selected: Press the corresponding function key  3, 4, 5, 6.

- Set the dosage for each HF output, using function keys  and .


### Note

An additional bar icon appears in the display showing where the selected dosage lies within the allowable range.

## Selecting activation mode foot pedal/hand switch/auto-COAG

An activation mode must be set for each output that is activated: with the foot pedal, the hand switch or automatically.


Icons on menu bar activation mode:

 on	 off	Description
		HF activation and switch-off via the foot pedal
	auto COAG	HF activation via the foot pedal and automatic HF switch-off
		HF activation and switch-off via the hand switch
	auto COAG	HF activation via the hand switch and automatic HF switch-off
	auto COAG	Automatic switch-on/off for coagulation (only available in bipolar mode); Note - only available in bipolar mode - the maximum dosage setting is reduced to 50

- Press the function key .

The menu bar activation mode is displayed:


- Press function key  repeatedly until the desired activation mode is highlighted in black.
- Confirm the setting by pressing function key  and return to the main menu.
  - or -
- To exit the menu without adopting the changes: Press the function key  or wait for 30 seconds.

#### Note

If a monopolar output is allocated with an activation mode with auto-COAG off that is not supported by the selected operational mode, the system automatically selects the operational mode CONTACTCOAG.


If the activation mode foot pedal is selected, another menu bar is displayed:


Three settings can be chosen in this menu bar:

- activation of CUT and COAG
- activation of CUT only
- activation of COAG only

#### Option: Double foot pedal operates two different HF outputs

- To control two different HF outputs with a double foot pedal: Press function key  at both outputs repeatedly until the required function is displayed.
- Confirm the setting by pressing function key  and return to the main menu.
  - or -
- To exit the menu without adopting the changes: Press the function key  or wait for 30 seconds.

Display color for CUT: yellow

Display color for COAG: blue


#### Danger for patient and user from accidental activation of an output!

- Exercise extreme caution when using a double foot pedal.
- Listen for the acoustic signal for CUT or COAG.
- Avoid contact between tissues and the deactivated instrument.

#### Automatic disabling of the activation mode

If an activation mode that has already been selected (foot pedal or autoCOAG) is allocated to an output, the device will react as follows:

- New allocation of an activation mode is transferred to the output to be operated
- The icon of the activation mode on the card belonging to the original allocation is marked with a red cross, meaning the activation mode is disabled for the non-operational output.

## Special settings (Menu Special Functions)

### Instrument detection settings

The instrument detection system recognizes Aesculap accessories that carry a code (e.g. argon handle GN 380).

#### Note

The *Nelson<sup>deluxe</sup>* (GN 640) will confirm by an acoustic signal that an instrument with instrument detection was plugged in and that the instrument data are being transferred.

- Press function key  2.

The display shows:


#### Note

After a waiting time of 30 seconds, or when the user presses function key , the system returns to the main menu.

- Press function key  5.

The menu Special Functions is displayed:

- Set instrument recognition , using function key , as shown on the display.
- Confirm the setting by pressing function key  and return to the main menu.
- or -
- To exit the menu without adopting the changes:  
Press the function key  or wait for 30 seconds.

The following settings are available:

- Off: Instrument recognition deactivated
- Default setting: The HF device calls up the factory settings stored in the accessory. These settings can be changed manually, up to the maximum value for the respective accessory.
- Last setting: The HF device calls up the setting last adjusted for the accessory concerned.

#### Note

When an argon control unit GN 370 is connected to the *Nelson<sup>deluxe</sup>* (GN 640), the instrument recognition also affects the parameter settings at the argon control unit GN 370.

### autoCOAG activation delay settings

If the activation mode autoCOAG on/off is selected for a bipolar output, the autoCOAG activation delay can be used for setting the time interval between tissue contact and HF activation.

- Press function key  2.

The display shows:


#### Note

After a waiting time of 30 seconds, or when the user presses function key , the system returns to the main menu.

- Press function key  5.

The menu Special Functions is displayed:

- Set activation delay , using function key .

The autoCOAG activation delay can be adjusted in 10 steps:

Step	1	2	3	4	5	6	7	8	9	10
Time	0 s	0.2 s	0.4 s	0.6 s	0.8 s	1.0 s	1.5 s	2.0 s	3.0 s	4.0 s

As soon as step 10 is reached, another touch of the function key returns the setting to step 1.

- Confirm the setting by pressing function key  and return to the main menu.

- or -

- To exit the menu without adopting the changes:  
Press the function key  or wait for 30 seconds.

As soon as there is tissue contact, the *Nelson<sup>deluxe</sup>* (GN 640) sounds a short acoustic signal. When the HF output is activated at the end of the selected activation delay time, an acoustic signal sounds for the duration of the activation.

While the instrument is in contact with tissue, the background of the COAG settings range for the selected output is highlighted in blue.

#### Note

The autoCOAG activation delay allows the operating surgeon time for preparation without the HF active, for the duration of the delay time.

### Automatic deactivation settings

When automatic deactivation is selected, the device detects the end of coagulation and deactivates high frequency mode.

- Press function key  2.

The display shows:



#### Note

After a waiting time of 30 seconds, or when the user presses function key , the system returns to the main menu.

- Press function key  5.

The menu Special Functions is displayed:

- Set automatic deactivation , using function key , as shown on the display.
- Confirm the setting by pressing function key  and return to the main menu.  
- or -
- To exit the menu without adopting the changes:  
Press the function key  or wait for 30 seconds.

The coagulation gradient for deactivation can be adjusted in two steps and applies equally to all monopolar and bipolar outputs:

- Long: A higher coagulation gradient is reached until deactivation (forced coagulation)
- Short: A lower coagulation gradient is reached until deactivation (soft coagulation)

If the function autoCOAG on/off is selected, the maximum dosage is limited to 50. At the Long setting, the deactivation behavior changes for the bipolar outputs. The HF power is activated as soon as the forceps touches tissue; it is deactivated as soon as the forceps does not touch any tissue anymore.

#### Note

The higher the dosage setting, the quicker the deactivation threshold will be reached, and the more superficial is the coagulation.

WARNING

The setting autoCOAG on/off should not be used in endoscopic applications, since there would be a risk of accidental activation, e.g. when the instrument is introduced through the trocar.

### Selecting the argon function

The argon control unit is connected to the HF device via the CAN bus interface 18. Activation signals at the handle or foot pedal are transmitted to the argon control unit.

- Press function key  2.

The display shows:


#### Note

After a waiting time of 30 seconds, or when the user presses function key , the system returns to the main menu.

- Press function key  5.

The menu Special Functions is displayed:

- To select the argon function , press function key  as shown on the display.

Communication between the HF device and the argon control unit:

Communication correct: 

Communication interrupted: 

- Confirm the setting by pressing function key  and return to the main menu.  
- or -
- To exit the menu without adopting the changes:  
Press the function key  or wait for 30 seconds.

A monopolar output is assigned to the argon function. The following settings are available:

- OFF: The argon function is deactivated
- Monopolar 1 (argon function at output 9)
- Monopolar 2 (argon function at output 10)

## Memory reset

The seven separately programmable memories 1–7 can be reset to default settings.

- Press function key  2.

The display shows:


### Note

After a waiting time of 30 seconds, or when the user presses function key , the system returns to the main menu.

- Press function key  5.

The menu Special Functions is displayed:

- Memory reset:  
With the function key  call up the Restore default memory settings.
- Confirm the setting by pressing function key  and return to the main menu.  
- or -
- To exit the menu without adopting the changes:  
Press the function key  or wait for 30 seconds.

## Selecting a language

The device can be set to one of five different languages:

- German
- English
- French
- Spanish
- Italian

- Press function key  2.

The display shows:


### Note

After a waiting time of 30 seconds, or when the user presses function key , the system returns to the main menu.

- Press function key  5.

The menu Special Functions is displayed:

- Set language, using function key .
- Confirm the setting by pressing function key  and return to the main menu.  
- or -
- To exit the menu without adopting the changes:  
Press the function key  or wait for 30 seconds.

## Saving device settings


The HF device features seven memories, to which device settings (operational modes, activation modes, dosage settings etc.) can be saved individually.

- Change to the menu Memory by pressing function key  1. Seven individually allocatable memories are available there.

- To call up one of the seven memories: Briefly (< 2 s) press the corresponding function key  14.
- To save the present device settings in one of the seven memories: Hold down the corresponding function key  14 for at least 2 seconds.

The menu Name Setting is displayed:

- Use the function keys for entering the memory name (available characters: A–Z, 1–9, special characters, blank spaces).
- Confirm the name of the memory with the function key  and return to the main menu.

The name and number of the memory  appear at the upper left side of the display. If at this point the saved settings are altered, the memory name and number will be deleted.

- To interrupt the Save function: Press the function key  or wait for 30 seconds. The previous data remains in the memory.

## Calling up suggested (default) dosages

Device settings concerning certain applications are stored in another seven memories. Starting from these settings, the parameters can be changed individually. However, the content of the preset memories will not be changed.

Fixed default settings of the seven preset memories:

- GEN-SURG – Memory 8: General surgery
- GYN – Memory 9: Gynecology (Activation is divided between two outputs!)
- NEURO – Memory 10: Neurosurgery
- ARTHRO – Memory 11: Arthroscopy
- UROLOGY – Memory 12: Urology
- ARGON-O – Memory 13: Argon function for General Surgery
- ARGON-E – Memory 14: Argon function for endoscopy

**Risk for patients and users due to accidental activation of an output when using memory GYN with a double foot pedal!**

- Exercise extreme caution when using a double foot pedal.
- Listen for the acoustic signal for CUT or COAG.
- Avoid contact between tissues and the deactivated instrument.

- Change to the menu Memory by pressing function key  1.
- To display an overview of the seven preset memories: Press function key  2.

- To call up one of the seven memories: Press the corresponding function key  14.

The name and number of the memory  appear at the upper left side of the display. If at this point the device settings are altered, the memory name and number will be deleted.

## Connecting the accessories

Combinations of accessories that are not mentioned in the instructions for use may only be employed if they are appropriate for the intended application and if they do not negatively impact performance parameters or safety requirements.

The insulation on the accessories (e.g., HF cables, instruments) must be sufficient for their maximum peak output voltage (see IEC 60601-2-18).

- Do not connect the accessory cable before the device has been switched on.
- Connect the foot pedal plug at foot pedal connection socket 11.
- Connect the plug of the neutral electrode cable at neutral electrode connection socket 12.
- Connect the potential equalization to potential equalization connector 22.
- To assign the data interface 18: Only connect cables and accessories identified by Aesculap as intended to be used for this purpose. (For specifications, see the instructions for use of the accessory, e.g. argon control unit GN 370).

Please address your Aesculap partner or Aesculap Technical Service with any inquiries in this respect.

### Activating the high frequency

- Be sure that the patient is prepared in such a way that the HF device can be applied without any risk to the patient.
- Verify that the accessories required are connected correctly.
- Be sure that all settings on the HF device are in conformance with the application requirements.
- Activate the required HF output via the hand switch, the foot pedal or through autoCOAG.

For each activation mode, a different acoustic signal sounds and the color of the corresponding field in the display changes:

CUT: Signal for CUT, display color yellow

COAG: Signal for COAG, display color blue

To warn of the approaching maximum time of HF activity, the volume of the acoustic signal rises to 100 % 7 seconds before the HF current is switched off automatically.


#### Note

If a Vario pedal (GN 323 or GN 325) is used for coagulation, a bar symbol appears on the display  indicating how far the Vario pedal is pressed down and what percentage of preprogrammed power has been released.

### 4.3 Functional test

Prior to each use, carry out a functional test of the device:

- Be sure that the accessories used for the functional test are not damaged in any visually detectable way.
- Prepare and set up the device, see Preparation and setup.
- Check the functionality of the following elements, one after the other, in the following sequence:
  - Power ON switch , indicator lamp (ON)
  - Automatic self-test after switching on: Acoustic signal ("sound check") and display ("display test")

- Function keys
- Correct operation (see Operating the device)
- Memory
- Bipolar activation
- Monopolar activation
- Neutral electrode monitoring (error message on display if neutral electrode not plugged in correctly)
- Mains power OFF switch 

## 5. Cleaning, disinfecting, sterilization and product care

### Note

Accessories such as the electrode handle, active electrodes, neutral electrodes and cables must be prepared, applied and sterilized as described in their corresponding instructions for use.

*Do not sterilize single-use electrodes, disposable cables and single-use handles!*

### 5.1 Disinfecting and cleaning

#### Risk of electric shock and burns!

- Be sure to unplug the device before cleaning it.
  - Do not use flammable or explosive cleaning or disinfecting solutions.
  - Be sure that no fluids are allowed to penetrate the device.
- Only use cleaning and disinfectant solutions that are allowable for cleaning surfaces under national laws.
  - Never lay the device in water or cleaning solution.
  - Never disinfect the product mechanically or by boiling it.
  - Allow any fluid that penetrates the product to run out immediately.
  - Apply cleaning and disinfectant solution.
  - Wipe away cleaning and disinfectant solution with a sponge or cloth that has been moistened in clean water.
  - Dry the unit with a clean, lint-free cloth.

### 5.2 Care and checking for damage

- Inspect all cables, and in particular electrode cables, for any damage to their insulation.
- Perform a safety inspection annually (see Safety inspection).

### 5.3 Sterilization method and parameters

The high-frequency surgical device **Nelson**<sup>deluxe</sup> (GN 640) must not be sterilized.

## 6. Maintenance

The high-frequency surgical device **Nelson**<sup>deluxe</sup> (GN 640) is a maintenance-free product.

### 6.1 Safety inspection

A safety inspection must be carried out annually.

- Only have the product and its accessories inspected by persons who possess the requisite training, knowledge and experience and who are authorized to do so.

The person performing the inspection documents the results of the inspection as well as the measurement values in conformance with the inspection protocol form.

- In the event significant deviations from the inspection protocol values are found, or if the prescribed maximum values are exceeded: Send the device to the service facility indicated herein (see Technical Service).

AESCULAP <sup>®</sup>					Inspection protocol – safety inspection		
INSPECTION PERFORMED:					annually		
TYPE OF DEVICE:					High-frequency surgical device (GN 640) from serial number 2000		
MANUFACTURED BY:					AESCULAP AG & CO. KG, Am Aesculap-Platz, 78532 Tuttlingen/Germany		
SN: .....					INVENTORY NO.: .....		
OPERATOR: .....							
SCOPE OF INSPECTION					Finding		
1. VISUAL INSPECTION					Satisfactory	Unsatisfactory	
1.1 Inspection of power cable for visually detectable damage							
1.2 Inspection of foot pedal for visually detectable damage							
1.3 Inspection of coagulation cable for visible signs of damage (record their number); check the fitting of the plug connector							
1.4 Rating plate still present and readable							
1.5 Test to ensure device achieves default values							
1.6 Overall condition of device: Dirt, damage							
2. ELECTRICAL INSPECTION					Measured value	OK Yes/No	
2.1 LF ground leakage current s.f.c. according to DIN EN 60601-1					$I_{\max} = 1 \text{ mA}$		
2.2 LF substitute device leakage current acc. to VDE 0751-1					$I_{\max} = 1 \text{ mA}$		
2.3 LF patient leakage current s.f.c. acc. to DIN EN 60601-1					$I_{\max} = 0.05 \text{ mA}$		
2.4 LF substitute patient leakage current with mains voltage at the application unit acc. to VDE 0751-1					$I_{\max} = 0.05 \text{ mA}$		
2.5 HF leakage current monopolar DIN EN 60601-2-2					$I_{\max} = 150 \text{ mA}$		
2.6 HF leakage current bipolar DIN EN 60601-2-2					$I_{\max} = 71 \text{ mA}$		
2.7 Insulation resistance (test voltage = 500 V DC)							
2.7.1 Mains against HF output					$R_{\min} = 7 \text{ M}\Omega$		
2.7.2 Mains against housing					$R_{\min} = 2 \text{ M}\Omega$		
2.7.3 HF output against housing							
1. Bipolar HF output against housing					$R_{\min} = 2 \text{ M}\Omega$		
2. Bipolar HF output against housing					$R_{\min} = 2 \text{ M}\Omega$		
1. Monopolar HF output against housing					$R_{\min} = 2 \text{ M}\Omega$		
2. Monopolar HF output against housing					$R_{\min} = 2 \text{ M}\Omega$		
2.8 Ground wire resistance incl. power cord acc. to VDE 0751-1					$0,3 \Omega$ at $\geq 0.2 \text{ A}$		
2.9 HF output with induction-free resistance							
	HF output	Operational mode	Dos- age	$R_{\text{load}}$	Rated value		
Monopolar		MICROCUT 1/2	300	500 [ $\Omega$ ]	300 watt $\pm$ 20 %		
		POWERCUT 1	300	500 [ $\Omega$ ]	300 watt $\pm$ 20 %		
		POWERCUT 2	250	500 [ $\Omega$ ]	250 watt $\pm$ 20 %		
		POWERCUT 3	200	500 [ $\Omega$ ]	200 watt $\pm$ 20 %		
		POWERCUT 4	150	500 [ $\Omega$ ]	150 watt $\pm$ 20 %		
		SEALCUT	100	150 [ $\Omega$ ]	100 watt $\pm$ 20 %		
		CONTACTOCOAG	120	500 [ $\Omega$ ]	120 watt $\pm$ 20 %		
		SOFTCOAG	100	300 [ $\Omega$ ]	100 watt $\pm$ 20 %		
		SPRAYCOAG	120	800 [ $\Omega$ ]	120 watt $\pm$ 20 %		
		FORCEDCOAG	120	1 000 [ $\Omega$ ]	120 watt $\pm$ 20 %		
Bipolar		BICUT 1	100	600 [ $\Omega$ ]	100 watt $\pm$ 20 %		
		BICUT 2	100	100 [ $\Omega$ ]	100 watt $\pm$ 20 %		
		BICOAG	100	100 [ $\Omega$ ]	100 watt $\pm$ 20 %		
2.10 DC resistance between the two HF output poles, DIN EN 60601-2-2; section 59.105							
2.11 Function test of the monitoring circuit for the neutral electrode, DIN EN 60601-2-2, section 59.101					Performed		
2.12 Function test according to instructions for use					Performed		
Place/Date		Inspector		Operator			

## 7. Recognizing and remedying errors

### 7.1 Calling up the activation, error and On Phase memories

- Press function key  2.

The display shows:


Note

After a waiting time of 30 seconds, or when the user presses function key , the system returns to the main menu.

- Press function key  4 to call up the menu History Record.

A pop-up menu is displayed.


- Press function key  to select the required memory menu.

- Confirm the setting by pressing function key .

- or -

- Leave the menu:

Press the function key  or wait for 30 seconds.

- Keys  are used for going from page to page within the menu; keys  are pressed to jump to the beginning or the end of the list.

#### Activation memory

The activation memory saves the last 600 activations of the HF device. The following parameters are shown:

- Date of activation (DD.MM.YY)
- Time of activation (hh:min)
- Duration of activation (in seconds)
- Output
  - BI1: bipolar 1
  - BI2: bipolar 2
  - MO1: monopolar 1
  - MO2: monopolar 2
- Dosage

- Mode
  - CUT0: POWERCUT 1/BICUT 1 COAG0: CONTACTCOAG / BICOAG
  - CUT1: POWERCUT 2/BICUT 2 COAG1: SOFTCOAG
  - CUT2: POWERCUT 3 COAG2: SPRAYCOAG
  - CUT3: POWERCUT 4 COAG3: FORCEDCOAG
  - CUT4: MICRO CUT 1
  - CUT5: MICRO CUT 2
  - CUT6: SEALCUT
- Activation mode on/off (H: hand, A: automatic, F: foot pedal)
- I<sup>2</sup>t-value (in 0.1 A<sup>2</sup>s)
- Instrument used
  - unknown: Instrument without instrument recognition
  - seven-digit number: 3rd to the 5th digit of the article number, and the last 4 digits of the serial number of the instrument

Activation							History Record		
date	time	d	out	dose	mode	on-off	I <sup>2</sup> t	instr	
21.08.03	11:03:37	02	BI1	9,0	COAG0	F-F	000	unknown	
21.08.03	11:03:34	02	BI1	40	CUT1	F-F	000	unknown	
21.08.03	11:03:24	04	BI2	40	CUT0	F-F	000	unknown	
21.08.03	11:03:19	04	BI2	40	COAG0	F-F	000	unknown	
21.08.03	11:03:10	03	MO1	40	COAG3	F-F	000	unknown	
21.08.03	11:03:05	03	MO1	80	CUT0	F-F	000	unknown	
21.08.03	11:03:00	03	MO2	80	COAG2	H-H	000	unknown	
21.08.03	11:02:55	04	MO2	150	CUT0	H-H	000	unknown	
21.08.03	11:02:47	<1	MO2	80	COAG2	H-H	000	unknown	
21.08.03	11:02:12	<1	MO2	150	CUT0	H-H	000	unknown	

#### Error memory (Error by date; Error by number)

The error memory saves the last 100 error messages generated by the HF device. High-priority errors are saved for a longer period. The following parameters are shown:

- Date (DD.MM.YY)
- Time
- Error number
- Error text

Error by date				History Record	
date	time	No.	description		
21.08.03	11:17:47	11	NE disconnected (dynamic)		
21.08.03	11:17:43	10	NE short circuit		
21.08.03	11:17:33	12	NE disconnected (static)		
21.08.03	11:17:27	12	NE disconnected (static)		
21.08.03	11:17:21	8	Keyboard time exceeded		
21.08.03	11:15:52	48	PIN-Error monopolar		
21.08.03	11:14:49	56	Time exceeded		
21.08.03	11:14:09	9	Keyboard multiple key pressed		
21.08.03	11:01:52	47	Double key		
21.08.03	10:42:01	169	Hardware incompatible component		


## On Phase memory

The On Phase memory saves the last 100 on phases of the HF device. The following parameters are shown:

- Date (DD.MM.YY)
- Switch-on time
- Duration of on phase (hh:mm)
- Type of neutral electrode (1 = type 1, 2 = type 2, 3 = type 3)
- Maximum HF switching time
- Instrument recognition setting (0 = off, 1 = default, 2 = last settings)
- autoCOAG automatic switch-off setting (S = short, L = long)
- Argon allocation (0 = none, 1 = monopolar 1, 2 = monopolar 2)
- Highest occurred temperature of the monopolar power amplifier

On Phase		History Record						
start date	time	on time hh:mm				auto coag		
21.08.03	11:20:33	00:04	1	30	1	L	0	37
21.08.03	11:10:28	00:03	1	30	1	S	0	38
21.08.03	10:41:03	00:40	1	30	1	L	0	38
21.08.03	10:33:21	00:08	1	30	1	L	0	36
21.08.03	09:27:23	01:06	1	30	1	L	0	35
21.08.03	08:05:09	01:21	1	30	1	L	0	35
21.08.03	07:55:09	00:06	1	30	1	S	0	38
21.08.03	07:45:09	00:03	1	30	1	S	0	36
21.08.03	07:34:08	00:05	1	30	1	S	0	35
21.08.03	07:25:08	00:02	1	30	1	S	0	35

- Press function key  3 to call up the menu Service Functions.

- Press function key  to call up the required service function or to activate (ON) or deactivate (OFF) the service mode.
- Confirm the setting by pressing function key  and return to the main menu.
- or -
- To exit the menu without adopting the changes: Press the function key  or wait for 30 seconds.

## 7.2 Calling up the menu Service Functions

The following data can be retrieved in the Service Functions menu:

- Circuit board identification: shows art. no., serial no. and the status of any alterations made in the integrated circuit board
- Actual value: shows various analog and digital actual values in the microcontroller periphery (see Service Manual)
- Service mode: OFF: Error list (see Error messages on the display) ON: Error list (see Service Manual)

- Press function key  2. The display shows:


### Note







After a waiting time of 30 seconds, or when the user presses function key , the system returns to the main menu.

## 7.3 Error messages on the display







If the device detects errors in the course of a self-test, these errors are displayed with an error number, a pictogram and an error text.

Depending on the error detected, the device disables the monopolar or bipolar HF outputs or the monopolar and bipolar HF outputs.

No.	Pictogram	Description	Cause	Remedy
2		Switch-on Test handle activated	Handle key monopolar pressed when starting	Release handle key or disconnect handle cable
3		Switch-on Test handle activated	Handle key bipolar pressed when starting	Release handle key or disconnect handle cable
4		Switch-on Test pedal activated	Foot pedal key pressed when starting	Release or disconnect foot pedal
5		Switch-on Test keyboard activated	Plastic foil keyboard pressed when starting	Release key
8		Keyboard time exceeded	Plastic foil keyboard held pressed for too long	Release key Replace plastic foil keyboard
9		Keyboard multiple key pressed	Multiple keys pressed	Release key Replace plastic foil keyboard
10		NE short circuit	Neutral electrode: Error from short circuit	Check NE cable Check NE cable plug-in connection
11		NE disconnected (dynamic)	Neutral electrode triggered (dynamic)	Check entire surface of NE application If okay, generate new reference values by disconnecting the plug and then connecting it again
12		NE disconnected (static)	Neutral electrode triggered (static) or connection severed	Check NE cable Check NE cable plug-in connection
20		Pedal information	Foot pedal position recognized incorrectly – monopolar	Check foot pedal
21		Pedal type	Foot pedal type recognized incorrectly – monopolar	Check foot pedal
22		Pedal type <-> selection	Foot pedal type recognized different from default pedal type – monopolar	Change pedal type setting on the device to match that of the foot pedal being used
23		Pedal type <-> AutoCOAG	Vario function of the Vario pedals (GN 323, GN 325) is not supported in autoCOAG mode – monopolar	Error message is issued only at the first actuation of the pedal. After that, the Vario pedals (GN 323, GN 325) behave like the dual foot pedals (GK 223, GK 324), i. e. the HF output operates at 100 % power both in CUT mode and in COAG mode
30		Pedal information	Foot pedal position recognized incorrectly – bipolar	Check foot pedal
31		Pedal type	Foot pedal type recognized incorrectly – bipolar	Check foot pedal
32		Pedal type <-> selection	Foot pedal type recognized different from default pedal type – bipolar	Change pedal type setting on the device to match that of the foot pedal being used

No.	Pictogram	Description	Cause	Remedy
33		Pedal type <-> AutoCOAG	Vario function of the Vario pedals (GN 323, GN 325) is not supported in autoCOAG mode – bipolar	Error message is issued only at the first actuation of the pedal. After that, the Vario pedals (GN 323, GN 325) behave like the dual pedals (GK 223, GK 324), i. e. the HF output operates at 100 % power both in CUT mode and in COAG mode
41		Incorrect output	Dosage monitoring (voltage) – monopolar	Switch device off and then start it again
42		Incorrect output	Dosage monitoring (current) – monopolar	Switch device off and then start it again
43		Incorrect output	Dosage monitoring (HF-on) – monopolar	Switch device off and then start it again
44		AutoCOAG	Criteria for automatic HF deactivation not achieved (time exceeded)	Switch to foot pedal activation
45		Temperature monitoring	Allowable temperature exceeded at monopolar end stage	Allow device to cool down
46		Time exceeded	Maximum allowable HF activation time exceeded – monopolar	Release foot pedal/handle
47		Double key	Two keys pressed on handle – monopolar	Release keys Check handle for possible short circuit
48		PIN-Error monopolar	Plug put in false socket – monopolar	Use correct socket (see chapter 4.1)
51		Incorrect output	Dosage monitoring (voltage) – bipolar	Switch device off and then start it again
52		Incorrect output	Dosage monitoring (current) – bipolar	Switch device off and then start it again
53		Incorrect output	Dosage monitoring (HF-on) – bipolar	Switch device off and then start it again
54		AutoCOAG	Criteria for automatic HF deactivation not achieved (time exceeded)	Switch to foot pedal activation
55	-	Output impedance	Incorrect output impedance – bipolar	Check HF cable and instrument Only use recommended accessories
56		Time exceeded	Maximum allowable HF activation duration exceeded – bipolar	Release foot pedal/handle
57		Double key	Two keys pressed on handle – bipolar	Release keys Check handle for possible short circuit
60		Hardware – CAN	Faulty communication – CAN	Check CAN interface cable

No.	Pictogram	Description	Cause	Remedy
61		Hardware – RS 232	Faulty communication – RS 232	Check RS 232 interface cable
80	-	External unit	Error message generated by external device at CAN-Bus and RS 232 interface	Check for error message on external device
90		Risk of burns	Risk of burns at the neutral electrode owing to high output (I <sup>2</sup> *t monitoring)	Check neutral electrode application area, lower dosage
95		Temperature monitoring	Excessive temperature in power supply	Allow device to cool down
100		Hardware	Hardware error	Switch device off and then start it again
102		Hardware-watch	Error in time base comparison	Switch device off and then start it again Change battery (see Service Manual)
107		Instrument detection	Error in bipolar instrument detection	Check HF cable and instrument Only use recommended accessories
108		Instrument detection	Error in monopolar instrument detection	Check HF cable and instrument Only use recommended accessories
109		Instrument detection	APC probe not connected	Connect APC probe
116		Hardware protective circuit bipolar	Bipolar instrument and/or cable defective	Replace bipolar instrument and/or cable
117		Instrument detection – wrong values	Instrument detection bipolar – data incompatible	Use only instruments with instrument detection, suited for the <i>Nelson</i> <sup>deluxe</sup> (GN 640) beginning with SN 2000
118		Instrument detection – wrong values	Instrument detection monopolar – data incompatible	Use only instruments with instrument detection, suited for the <i>Nelson</i> <sup>deluxe</sup> (GN 640) beginning with SN 2000
150		Battery	Empty battery	Change battery
200		Software	Software	Switch device off and then start it again

## 7.4 Changing fuses

Allowable fuses:

2 pcs. Aesculap art. no. TA 021 473:

G melting fuse, time-delay (T) 6.3 A,  
breaking capacity (1 500 A with 250 V/50–60 Hz)

- Pull out the device plug before changing the fuse!
- Use a small screwdriver to release the clips at the fuse holder **20**.
- Pull out the fuse holder.
- Change both fuses.
- Reinsert the fuse holder in such a way that it audibly snaps into place.

*Note*

*If fuses burn out frequently, the device is malfunctioning and must be repaired (see Technical Service).*

## 8. Technical Service

For service, maintenance or repairs contact your Aesculap/B. Braun representative in your country.

Carrying out modifications on medical equipment renders the guarantee, and any licenses, null and void.

### Service addresses

AESCULAP Technischer Service  
Am Aesculap-Platz  
78532 Tuttlingen/Germany  
Phone: +49 (7461) 95 27 00  
Fax: +49 (7461) 16 28 87  
E-mail: ats@aesculap.de

or in the US:

AESCULAP Inc.  
Attn. AESCULAP Technical Services  
615 Lambert Pointe Drive  
Hazelwood  
MO, 63042  
AESCULAP Repair Hotline  
Phone: +1 (800) 214-3392  
Fax: +1 (314) 895-4420

Other service addresses can be obtained from the address indicated above.

## 9. Disposal

*Note*

*The user institution is obliged to process the product before its disposal, see Cleaning, disinfecting, sterilization and product care.*

*Always adhere to national regulations when disposing of or recycling the product or its components!*

Products carrying this symbol are subject to separate collection of electric and electronic devices. Within the European Union, disposal is taken care of by the manufacturer as a free-of-charge service.

In case of any questions concerning disposal of the product, please contact your national B. Braun/Aesculap agency, see Technical Service.

## 10. Accessories and spare parts

Designation	Art. no.
Fuse	TA 021 473
Material for stacking devices <b>without</b> integrated stacking equipment (e.g., GN 370)	GN 370 830
Material for stacking devices <b>with</b> integrated stacking equipment (e.g., GN 090)	GN 370 831
Undivided neutral electrode (type 1)	e.g. GK 105
Divided disposable neutral electrode (type 2)	e.g. GK 106
Reusable divided neutral electrode (type 3)	GN 320

Information about additional accessories and replacement parts can be found in Aesculap brochure C-304-81.

## 11. Technical specifications

Model	<i>Nelson</i> <sup>®</sup> deluxe (GN 640)	
Voltage range	100–240 V	
Power consumption	4.7 A with 100 V 3.6 A with 120 V 1.9 A with 220 V 1.7 A with 240 V	
Protection class according to EN 60601-1	I	
Classification according to EG Directive 93/42/EEC	IIb	
Foot pedal circuit	fused in conformance with IEC 60601, certified for use in a "medical environment"	
HF output – monopolar	MICROCUT 1 and 2	300 W at 500 Ohms
	POWERCUT 1	300 W at 500 Ohms
	POWERCUT 2	250 W at 500 Ohms
	POWERCUT 3	200 W at 500 Ohms
	POWERCUT 4	150 W at 500 Ohms
	SEALCUT	100 W at 150 Ohms
	SOFTCOAG	100 W at 300 Ohms
	CONTACTCOAG	120 W at 500 Ohms
	FORCEDCOAG	120 W at 1 000 Ohms
	SPRAYCOAG	120 W at 800 Ohms
HF output – bipolar	BICUT 1	100 W at 600 Ohms
	BICUT 2	100 W at 100 Ohms
	COAG	100 W at 100 Ohms
Frequency	447 kHz	
Modulation frequency – monopolar	20 kHz	
Modulation frequency – bipolar	1 kHz	
Operational mode	Int 10 s/30 s	
Miniature fuse	T 6.3 A Breaking capacity: 1 500 A at 250 V/50–60 Hz	
Weight	8.6 kg	
Dimensions (W x H x D)	305 x 175 x 305 mm	
Type of application part according to EN 60601-1	CF	
EMV	IEC 60601-1-2	
Standards	IEC 60601-1 IEC 60601-2-2	

Standards	CSA-CUS  LR 50 151
	

## 12. Index

### A

- Accessories 33, 40, 47, 51, 52
  - Connecting the accessories 50
- Activation
  - Accidental activation 34, 46, 48, 50
  - Activating high frequency mode 51
- Activation memory 54
- Activation mode 43, 45, 54
  - Activation mode autoCOAG 47
  - Activation mode foot pedal 46
  - AutoCOAG 45
  - Automatic disabling of the activation mode 46
  - Foot control 45
  - Foot pedal 45
  - Hand switch 45
- Active electrodes 34
  - Connecting active electrodes 42
- Argon 47, 55
  - Argon control unit GN 370) 50
  - Selecting the argon function 48
- AutoCOAG
  - Activation mode autoCOAG 45, 47
  - Activation/deactivation via autoCOAG 37, 48
  - Automatic deactivation 48
- autoCOAG 36, 45, 55
  - autoCOAG activation delay 47
  - autoCOAG Long 48
  - autoCOAG Short 42, 48
  - Automatic deactivation 48

### C

- Changing fuses 59
- Circuit board identification 55
- Coagulation gradient 48

### D

- Date 42, 54, 55
- Diagrams
  - Adaptation curves 38
  - Maximum peak output voltage 40
  - Output diagrams 38
- Disposal
  - Disposal 32, 59
- Dosage 54
  - Dosage range 45
  - Dosage setting 41, 43, 45

Maximum dosage 45  
Minimum dosage 45

## E

Error 37, 54, 55  
Error list 56, 57  
Error memory 54  
Error message 37  
Error messages 55  
Error number 55  
Error text 55

Error memory 54

## F

Foot control  
Activation mode foot control 45  
Foot pedal 32, 34, 35, 48, 51  
Activation mode foot pedal 45, 46  
Connecting the foot pedal 42, 50  
Connection socket foot pedal 32, 42  
Double foot pedal operates two outputs 46  
Foot pedal plug 42

## H

Hand switch 34, 51  
Activation mode hand switch 45

## I

Instrument recognition 55  
Instrument recognition settings 47

## M

Maximum HF switching time 42  
Maximum time of HF activity 51  
Memory 51  
Activation memory 54  
Calling up memories 50  
Data memory functions 37  
Error memory 54  
Function key memory menu 32  
Individually allocatable memories 50  
Memory name setting 50  
Memory reset 49  
On Phase memory 54, 55  
Preset memories 50  
Saving device settings 50  
Menu Service Functions 55  
Menu Special Functions 47, 48, 49  
Menu User Settings 42, 43  
Monitoring functions 37  
Neutral electrode monitoring 37, 51  
Self-monitoring 55

## N

Neutral electrode 32, 34, 35, 37, 41, 52  
Application of the neutral electrode 34  
Connecting the neutral electrode 41  
Connection socket neutral electrode 32, 41  
Neutral electrode fault 41  
Neutral electrode monitoring 37, 51  
Neutral electrode symbol 41  
Neutral electrode type 1 – 37, 42, 55, 59  
Neutral electrode type 2 – 37, 42, 55, 59  
Neutral electrode type 3 – 35, 37, 42, 55, 59  
Plug of the neutral electrode cable 41, 50

## O

On Phase memory 54, 55  
Operating environment 33  
Operational mode 35, 45  
BICOAG – 37, 43, 45, 60  
BICUT 1 – 37, 43, 45, 60  
BICUT 2 45  
BICUT 2 – 37, 43, 60  
Bipolar coagulation 37, 38, 39, 40, 45  
Bipolar incisions 37, 38, 39, 40, 43, 44, 45  
CONTACTCOAG 36, 43, 45, 60  
FORCEDCOAG 36, 43, 45, 60  
MICROCUT 1 – 36, 43, 45, 60  
MICROCUT 2 – 36, 43, 45, 60  
Monopolar coagulation 36, 38, 40, 44, 45  
Monopolar incisions 36, 38, 40, 43, 44, 45  
POWERCUT 1 – 36, 43, 45, 60  
POWERCUT 2 – 36, 43, 45, 60  
POWERCUT 3 – 36, 45, 60  
POWERCUT 4 – 36, 45, 60  
SEALCUT 36, 43, 45, 60  
Selecting the operational mode COAG 44  
Selecting the operational mode CUT 44  
Selection of an operational mode 41, 43  
SOFTCOAG 36, 43, 45, 60  
SPRAYCOAG 36, 43, 45, 60  
Using the operational modes 35  
Output socket  
Output socket HF – bipolar 1 – 32, 42  
Output socket HF – bipolar 2 – 32, 42  
Output socket HF – monopolar 1 – 32, 42  
Output socket HF – monopolar 2 – 32, 42

## S


Selecting a language 49  
Selecting and de-selecting HF outputs 43  
Self-test 37, 51  
Service addresses 59  
Service mode 55

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